



ASSYST bvba / A.S.O.W. bvba
Hellegatstraat 13a - 2590 Berlaar - Belgium
Tel: +32 495 50 61 14 / +32 496 83 70 27
Website: www.assyst.org / www.artsuppliesonweb.com
Email: ao@assyst.org / vera.opsommer@assyst.org

CRITTER CLAY

DESCRIPTION:

Critter Clay is the number one choice of taxidermists and professionals worldwide because of its highly detailed, fine texture and low shrinkage. Critter Clay dries to a hard mass. It is a lightweight, extra strong, ready to use, water based wet modeling clay. Can be painted with all paints.

PROPERTIES:

- ✓ Selfhardening, no ovens needed
- ✓ 1-component, no mixing required
- ✓ Water based
- ✓ Sense of plasticity
- ✓ Gets rock hard
- ✓ Lightweight when fully cured
- ✓ Extra strong
- ✓ 0-1% shrinkage
- ✓ Much detail with fine structure
- ✓ Safe, harmless
- ✓ Gluten-free
- ✓ Few or no fault lines
- ✓ Life time can be extended by freezing and defrosting the product
- ✓ Easy to use, paint and finish
- ✓ Easy to clean with water

INSTRUCTIONS FOR USE:

- ✓ Ready to go from the bag.
- ✓ Once exposed to the air, Critter Clay begins to harden.
- ✓ Use water to wet your fingers and tools to better shape and add detail.
- ✓ Air drying and curing of the sculpture varies with the size of the project and the conditions you are working in. (Drying and curing takes an average of 2-6 days depending on thickness, size and environment).
- ✓ Keep Critter Clay in a plastic bag as long as you don't want it to start to harden.
- ✓ For extra shelf life, you can freeze and thaw Critter Clay.
- ✓ Critter Clay accepts all paints.
- ✓ Safe for children to use under adult supervision.

VARIOUS APPLICATIONS:

- ✓ Taxidermy (set glass eye, set animal ear, ...)
- ✓ In the classroom and at camp (handprints of children, modeling, ...)
- ✓ Ceramic prototyping
- ✓ Making of pots/vases, ...
- ✓ Endless applications for artists, model makers, students, ...

TIPS:

- ✓ When the clay begins to dry slightly by leaving the bag open, you can rewet the Critter Clay, and leave it in the bag overnight. The less water you use, the better the results you will get.